

Målinger er ca. og vejledende til brug for sammenligning af de enkelte produkters effekt i forhold til behandlingstid/afstand.		10.000 lux Behandlingstid < = ½ time	5000 lux Behandlingstid Ca. ½ - ¾ time	2500 lux Behandlingstid Ca. 1 - 2 timer
1. Innosol RONDO	55 w			
4000 K 840		22,5 cm	33,5 cm	45 cm
2. Innosol LUCIA 72 W & Aurora 72 W	72 w			
4000 K 840		30 cm	45 cm	60 cm
5500 K 955		32 cm	47 cm	65 cm
3. Innosol MESA 110 & Lucia Mega 110 W & Origo 110 W	110 w			
4000 K 840		42 cm	63 cm	84 cm
5500 K 955		44 cm	65 cm	86 cm
6500 K 865		45 cm	66 cm	87 cm
4. Innosol MESA 160	160 w			
4000 K 840		52,5 cm	79 cm	105 cm
5500 K 955		54 cm	82 cm	107 cm
5. Innosol MESA Mega Bright	160 w			
4000 K 840		57,5 cm	86 cm	115 cm
5500 K 955		59 cm	89 cm	118 cm
6. Innosol MESA Mega Bright new	160 w			
4000 K 840		70 cm	105 cm	140 cm
5500 K 955		72 cm	107 cm	143 cm
7. Innosol AURORA	70 w			
4000 K 840		40 cm	60 cm	80 cm
8. Innosol Multi 255	110 w			
3000 K 830		35,5 cm	53 cm	71 cm
6000 K 960		37 cm	55 cm	74 cm
6500 K 865		38 cm	57 cm	76 cm
9. Innosol Multi 439	156 w			
3000 K 830		42,5 cm	64 cm	85 cm
6000 K 960		44 cm	66 cm	87 cm
6500 K 865		45 cm	67 cm	88 cm
10. Innosol Multi 454	216 w			
3000 K 830		47,5 cm	71 cm	95 cm
6000 K 960		49 cm	73 cm	97 cm
6500 K 865		51 cm	75 cm	99 cm
11. PER3	LED	15 cm	22,5 cm	30 cm
12. Boston Twin	2x 24	21 cm	31,5 cm	42 cm
13. Biolamp	144 w	35 cm	60 cm	91 cm
14. Luxo 20W	20 w	7 cm	13 cm	22 cm



Strandvangsvej 6 - 3060 Espergærde

Tlf. 49 13 41 32 - Fax 49 13 41 23 - e-mail: info@sundcentret.dk

Se dagslyslamperne på : www.sundcentret.dk